

Diese Preisliste wurde mit tagesaktuellen Preisen am 9.2.2010 um 12:44 Uhr aus dem Internet generiert.

## **Price list "01.01.2010"**


### **Business Area**

### **Medical Technology**

Headquarters:

HUPFER<sup>®</sup> Metallwerke GmbH & Co. KG  
Dieselstraße 20  
D-48653 Coesfeld  
Tel: +49 (0) 2541/805-0  
Fax: +49 (0) 2541/805-369  
info@hupfer.de  
<http://www.hupfer.de>

Product group	Page No.
Transport Trolleys->Closed Transport Trolleys->Articles	1
Transport Trolleys->Open Transport Trolleys->Articles	3
Transport Trolleys->Accessories for Closed Transport Trolleys->Articles	4
Transport Trolleys->Accessories for Open Transport Trolleys->Articles	6
Transport Trolleys->Transport Trolleys SW and MSW->Transport Trolleys, MSW, Standard version	7
Transport Trolleys->Transport Trolleys SW and MSW->Transport Trolleys, MSSW, Heavy-duty version	8
Transport Trolleys->Transport Trolleys SW and MSW->Serving Trolley SW - Standard version	8
Transport Trolleys->Accessories for Transport Trolleys SW and MSW->Articles	12
Transport Trolleys->Basket Trolleys, Platform Trolleys, Transport Trolleys, Single Ended Trolleys->Articles	14
General Conditions of Sale and Delivery / General Conditions	

**Transport Trolleys -> Closed Transport Trolleys -> Articles**Order No. article €7504643 CTWS-6 965/751/1465 Stainless Steel Container transport trolley standard-model on request

Dimensions W x D x H mm:	965 x 751 x 1465 mm
-----------------------------	---------------------

7504644 CTWS-9 1249/751/1465 Stainless Steel Container transport trolley standard-model on request

Dimensions W x D x H mm:	1249 x 751 x 1465 mm
-----------------------------	----------------------

7504641 CTW-6 965/754/1435 Stainless Steel Container transport trolley on request

Inside dimensions in mm:	740 x 640 x 1077 mm
Dimensions W x D x H mm:	965 x 754 x 1435 mm

7504642 CTW-9 1249/754/1435 Stainless Steel Container transport trolley on request

Inside dimensions in mm:	994 x 640 x 1077 mm
Dimensions W x D x H mm:	1249 x 754 x 1435 mm

7504647 STWS-6 965/751/1465 Stainless Steel Sterile goods transport trolley standard-model on request

Dimensions W x D x H mm:	965 x 751 x 1465 mm
-----------------------------	---------------------

7504648 STWS-9 1249/751/1465 Stainless Steel Sterile goods transport trolley standard-model on request

Dimensions W x D x H mm:	1249 x 751 x 1465 mm
--------------------------	----------------------

7504844 Closed transport cart 6 StE STWS-6, standard version with one fold down wire shelf on request

7504645 STW-6 965/754/1435 Stainless Steel Sterile goods transport trolley on request



Inside dimensions in mm:	740 x 640 x 1077 mm
Dimensions W x D x H mm:	965 x 754 x 1435 mm

7504845 Closed transport cart 9 StE STWS-9, standard version with one fold down wire shelf on request



7504646 STW-9 1249/754/1435 Stainless Steel Sterile goods transport trolley on request



Inside dimensions in mm:	994 x 640 x 1077 mm
Dimensions W x D x H mm:	1249 x 754 x 1435 mm

7504809 Transport trolley ATWS 965/754/1450 Stainless Steel standard-model on request



Dimensions W x D x H mm:	965 x 754 x 1450 mm
--------------------------	---------------------

7504808 Transport trolley ATWS 1249/754/1450 Stainless Steel standard-model on request









Dimensions W x D x H mm:	1249 x 754 x 1450 mm
--------------------------	----------------------

7504807 Transport trolley ATW 965/754/1450 Stainless Steel on request



7504806	Transport trolley ATW 1249/754/1450 Stainless Steel	Dimensions W x D x H mm:	965 x 754 x 1450 mm	on request
		Dimensions W x D x H mm:	1249 x 754 x 1450 mm	


**Transport Trolleys -> Open Transport Trolleys -> Articles**

Order No.	article		€
7504225	OCTW-6 825/675/1367 Stainless Steel Container transport trolley open	910,00 € / Pc	
		Dimensions W x D x H mm:	825 x 675 x 1367 mm
7504211	OCTW-9 1155/675/1367 Stainless Steel Container transport trolley open	1100,00 € / Pc	
		Dimensions W x D x H mm:	1155 x 675 x 1367 mm
7504228	OSTW-6 825/675/1367 Stainless Steel Sterile goods transport trolley open	750,00 € / Pc	
		Dimensions W x D x H mm:	825 x 675 x 1367 mm
7504212	OSTW-9 1175/675/1367 Stainless Steel Sterile goods transport trolley open	840,00 € / Pc	
		Dimensions W x D x H mm:	1175 x 675 x 1367 mm
0119022	Gastronorm Trolley RWG/20 B L-tray insert Stainless Steel	679,00 € / Pc	
		Specification:	High version for baking trays 600 x 400 mm
		Number of compartments:	1
		Insert dimensions:	Baking trays 600 x 400 mm and ISO-Norm container 600 x 400 mm

Insert mode:	Lengthways
Number of insert rails:	20
Capacity:	20 x 600 x 400
Type of insert:	L-rails with push-through protection on both sides
Clearance:	75 mm
Chassis:	4 swivel castors, 2 of them with total brake, Ø 125 mm
Load capacity:	100 kg
Net weight:	22 kg
External dimensions W x D x H in mm:	533 x 696 x 1817 mm

**Transport Trolleys -> Accessories for Closed Transport Trolleys -> Articles**

Order No.	article	€
-----------	---------	---

7504555	Insert frame set StU 296/595/60 Stainless Steel e-polished appr. 3/4 extension	90,00 € / Pc
---------	--	--------------



Dimensions W x D x H mm:	595 x 296 x 58 mm
--------------------------	-------------------

7504213	From high-quality Stainless Steel, approx. 350 mm drawer length, with 2 insert rails, for container and baskets according to DIN StU Insert rails 600/84/41 (1 pair)	39,00 € / Pair
---------	---	----------------



Dimensions W x D x H mm:	600 x 84 x 41 mm
--------------------------	------------------

7505009	From high-quality Stainless Steel 2 mm material thickness 1 Unit = 2 pcs Supporting surface 60 mm U-Profile Rails 1 pair	31,00 € / Pair
---------	--	----------------



Dimensions W x D x H mm:	36 x 600 x 41 mm
--------------------------	------------------

7504373	From high-quality Stainless Steel 2 mm material thickness 1 Unit = 2 pcs Supporting surface 14 mm Insert shelf 300/600/15 Stainless Steel for applying onto insert rails	65,00 € / Pc
---------	--	--------------



Dimensions W x D x H mm:	DIN StE 300 x 600 x 15 mm
--------------------------	---------------------------

0102972	Shelf solid 685/600 Stainless Steel	on request
---------	-------------------------------------	------------



0102682 Shelf solid 100/60 N5 Stainless Steel



0102002 Stainless Steel Wire Shelf 685 x 600 mm

on request



0101782 Shelf wire 100/60 N5 Stainless Steel

98,00 € / Pc



7501378 Shelf wire 772/590 for STWS-6 Stainless Steel

115,00 € / Pc



7501379 Shelf wire 1056/590 for STWS-9 Stainless Steel

130,00 € / Pc



7504232 Gallery 3-sided for STW / CTW 9

210,00 € / Pc



No self-assembly, will be assembled from manufacturer

7504233 Gallery, 4-sides for STW / CTW 9

245,00 € / Pc



No self-assembly, will be assembled from manufacturer

7500763 Cylinder locking for CTW / STW

105,00 € / Pc



No self-assembly, will be assembled from manufacturer

7504379 Card holder 117/50 Stainless Steel for CTW and STW

31,00 € / Pc



No self-assembly, will be assembled from manufacturer

7504380 Card holder DIN A5 Stainless Steel for CTW and STW

28,00 € / Pc



DIN A 5  
No self-assembly, will be assembled from manufacturer  
7504381 Card holder DIN A4 Stainless Steel for CTW and STW 35,00 € / Pc




DIN A4  
No self-assembly, will be assembled from manufacturer  
7500895 Paper Clamp, Stainless Steel for solid transport trolleys 38,00 € / Pc



Dimensions W x D x H mm:	90 x 18 x 70 mm
-----------------------------	-----------------

7504715 From high quality Stainless Steel, for CTW, CTWS, STW and STWS  
No self-assembly, will be assembled from manufacturer  
Transport locking with leaded bracket  
From high-quality Stainless Steel, for CTW and STW on request  
No self-assembly, will be assembled from manufacturer

### **Transport Trolleys -> Accessories for Open Transport Trolleys -> Articles**

Order No.	article	€
0102972	Shelf solid 685/600 Stainless Steel	on request
		
0102682	Shelf solid 100/60 N5 Stainless Steel	66,00 € / Pc
		
0102002	Stainless Steel Wire Shelf 685 x 600 mm	on request
		
0101782	Shelf wire 100/60 N5 Stainless Steel	98,00 € / Pc
		
7500762	Swivelled Stop bar for OCTW	122,00 € / Pc
		

7500362 For self-assembly  
Push Handle for OCTW and OSTW

105,00 € / Pc



Round tube, no self-assembly, will be assembled from manufacturer

7504357 Profile rail 594/40/80 Stainless Steel with 2 welded hooks







49,00 € / Pc



Dimensions W x D x H mm:	594 x 60 x 80 mm
-----------------------------	------------------

For lateral sticking to OCTW and OSTW version

### Transport Trolleys -> Transport Trolleys SW and MSW -> Transport Trolleys, MSW, Standard version

Order No.	article	€				
7504128	Transport trolley MSW 6x4/2 Stainless Steel standard version	315,00 € / Pc				
	 <table border="1"> <tr> <td>Shelf dimensions in mm:</td> <td>600 x 400 mm</td> </tr> <tr> <td>Total dimensions W x D x H in mm:</td> <td>700 x 500 x 895 mm</td> </tr> </table>	Shelf dimensions in mm:	600 x 400 mm	Total dimensions W x D x H in mm:	700 x 500 x 895 mm	
Shelf dimensions in mm:	600 x 400 mm					
Total dimensions W x D x H in mm:	700 x 500 x 895 mm					
7504129	Transport trolley MSW 6x4/3 Stainless Steel standard version	350,00 € / Pc				
	 <table border="1"> <tr> <td>Shelf dimensions in mm:</td> <td>600 x 400 mm</td> </tr> <tr> <td>Total dimensions W x D x H in mm:</td> <td>700 x 500 x 895 mm</td> </tr> </table>	Shelf dimensions in mm:	600 x 400 mm	Total dimensions W x D x H in mm:	700 x 500 x 895 mm	
Shelf dimensions in mm:	600 x 400 mm					
Total dimensions W x D x H in mm:	700 x 500 x 895 mm					
7504001	Transport trolley MSW 8x5/2 Stainless Steel standard version	270,00 € / Pc				
	 <table border="1"> <tr> <td>Shelf dimensions in mm:</td> <td>800 x 500 mm</td> </tr> <tr> <td>Total dimensions W x D x H in mm:</td> <td>900 x 600 x 940 mm</td> </tr> </table>	Shelf dimensions in mm:	800 x 500 mm	Total dimensions W x D x H in mm:	900 x 600 x 940 mm	
Shelf dimensions in mm:	800 x 500 mm					
Total dimensions W x D x H in mm:	900 x 600 x 940 mm					
7504002	Transport trolley MSW 10x6/2 Stainless Steel standard version	350,00 € / Pc				
	 <table border="1"> <tr> <td>Shelf dimensions in mm:</td> <td>1000 x 600 mm</td> </tr> <tr> <td>Total dimensions W x D x H in mm:</td> <td>1100 x 700 x 940 mm</td> </tr> </table>	Shelf dimensions in mm:	1000 x 600 mm	Total dimensions W x D x H in mm:	1100 x 700 x 940 mm	
Shelf dimensions in mm:	1000 x 600 mm					
Total dimensions W x D x H in mm:	1100 x 700 x 940 mm					
7504005	Transport trolley MSW 8x5/3 Stainless Steel standard version	310,00 € / Pc				
	 <table border="1"> <tr> <td>Shelf dimensions in mm:</td> <td>800 x 500 mm</td> </tr> <tr> <td>Total dimensions W x D x H in mm:</td> <td>900 x 600 x 940 mm</td> </tr> </table>	Shelf dimensions in mm:	800 x 500 mm	Total dimensions W x D x H in mm:	900 x 600 x 940 mm	
Shelf dimensions in mm:	800 x 500 mm					
Total dimensions W x D x H in mm:	900 x 600 x 940 mm					
7504006	Transport trolley MSW 10x6/3 Stainless Steel standard version	460,00 € / Pc				
						

Shelf dimensions in mm:	1000 x 600 mm
Total dimensions W x D x H in mm:	1100 x 700 x 940 mm

7504009 Transport trolley MSW 8x5/4 Stainless Steel standard version

520,00 € / Pc



Shelf dimensions in mm:	800 x 500 mm
Total dimensions W x D x H in mm:	900 x 600 x 1260 mm

### Transport Trolleys -> Transport Trolleys SW and MSW -> Transport Trolleys, MSSW, Heavy-duty version

Order No. article

€

7504107 Transport trolley MSSW 8x5/2 Stainless Steel heavy-duty version

450,00 € / Pc



Shelf dimensions in mm:	800 x 500 mm
Total dimensions W x D x H in mm:	900 x 600 x 1020 mm

7504110 Transport trolley MSSW 10x6/2 Stainless Steel heavy-duty version

500,00 € / Pc



Shelf dimensions in mm:	1000 x 600 mm
Total dimensions W x D x H in mm:	1100 x 700 x 1020 mm

7504113 Transport trolley MSSW 8x5/3 Stainless Steel heavy-duty version

550,00 € / Pc



Shelf dimensions in mm:	800 x 500 mm
Total dimensions W x D x H in mm:	900 x 600 x 1020 mm

7504114 Transport trolley MSSW 10x6/3 Stainless Steel heavy-duty version

585,00 € / Pc



Shelf dimensions in mm:	1000 x 600 mm
Total dimensions W x D x H in mm:	1100 x 700 x 1020 mm

### Transport Trolleys -> Transport Trolleys SW and MSW -> Serving Trolley SW - Standard version

Order No. article

€

0124002 General Purpose Trolley SW 6x4/2 Stainless Steel

295,00 € / Pc



Number of shelves:	2
Shelf dimensions in mm:	600 x 400 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	538 mm

Height of top shelf:	800 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	80 kg
Net weight:	13 kg
External dimensions W x D x H in mm:	700 x 500 x 894 mm

0124012 General Purpose Trolley SW 6x4/3 Stainless Steel

335,00 € / Pc



Number of shelves:	3
Shelf dimensions in mm:	600 x 400 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	251 mm
Height of top shelf:	800 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	120 kg
Net weight:	16 kg
External dimensions W x D x H in mm:	700 x 500 x 894 mm

0111932 General Purpose Trolley SW 8 x 5 / 2 Chrome Nickel Steel

240,00 € / Pc



Number of shelves:	2
Shelf dimensions in mm:	800 x 500 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	590 mm
Height of top shelf:	853 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	80 kg
Net weight:	15 kg
External dimensions W x D x H in mm:	900 x 600 x 947 mm

0111962 General Purpose Trolley SW 8 x 5 / 3 Chrome Nickel Steel

289,00 € / Pc



Number of shelves:	3
Shelf dimensions in mm:	800 x 500 mm
Load per unit area / shelf:	40 kg

Working height between shelves:	278 mm
Height of top shelf:	853 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	120 kg
Net weight:	20 kg
External dimensions W x D x H in mm:	900 x 600 x 947 mm

0111992 General Purpose Trolley SW 8x5/4 Stainless Steel

510,00 € / Pc



Number of shelves:	4
Shelf dimensions in mm:	800 x 500 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	278 mm
Height of top shelf:	1166 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	150 kg
Net weight:	25 kg
External dimensions W x D x H in mm:	900 x 600 x 1260 mm

0112022 General Purpose Trolley SW 8x5/5 Stainless Steel

590,00 € / Pc



Number of shelves:	5
Shelf dimensions in mm:	800 x 500 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	278 mm
Height of top shelf:	1479 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	150 kg
Net weight:	30 kg
External dimensions W x D x H in mm:	900 x 600 x 1574 mm

0112052 General Purpose Trolley SW 10x6/2 Stainless Steel

440,00 € / Pc



Number of shelves:	2
Shelf dimensions in mm:	1000 x 600 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	590 mm
Height of top shelf:	853 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	80 kg
Net weight:	19 kg
External dimensions W x D x H in mm:	1100 x 700 x 947 mm

0112082 General Purpose Trolley SW 10x6/3 Stainless Steel

530,00 € / Pc




Number of shelves:	3
Shelf dimensions in mm:	1000 x 600 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	278 mm
Height of top shelf:	853 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	120 kg
Net weight:	26 kg
External dimensions W x D x H in mm:	1100 x 700 x 947 mm

0112112 General Purpose Trolley SW 10x6/4 Stainless Steel




660,00 € / Pc



Number of shelves:	4
Shelf dimensions in mm:	1000 x 600 mm
Load per unit area / shelf:	40 kg
Working height between shelves:	278 mm
Height of top shelf:	1166 mm
Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm
Total load capacity:	150 kg
Net weight:	33 kg

0112142	General purpose trolley SW 10x6/5 Stainless Steel		External dimensions W x D x H in mm:	1100 x 700 x 1260 mm	836,00 € / Pc
			Number of shelves:	5	
			Shelf dimensions in mm:	1000 x 600 mm	
			Load per unit area / shelf:	40 kg	
			Working height between shelves:	278 mm	
			Height of top shelf:	1479 mm	
			Chassis:	4 swivel castors, 2 of them with brake, Ø 125 mm	
			Total load capacity:	150 kg	
			Net weight:	40 kg	
			External dimensions W x D x H in mm:	1100 x 700 x 1573 mm	

### Transport Trolleys -> Accessories for Transport Trolleys SW and MSW -> Articles

Order No.	article	€
0114506	Gallery 3-sided for SW/SSW-60 with shelf size 600/400 Stainless Steel 50 mm total height, 40 mm above shelf	110,00 € / Pc
		
	Stainless Steel wire Ø 8 mm Gallery height: 50 mm Gallery suitable for all shelves of a serving trolley	
0114502	For shelf-dimension: 600 x 400 mm Gallery 3-sided for SW/SSW-80 with shelf size 800/500 Stainless Steel 50 mm total height, 40 mm above shelf	115,00 € / Pc
		
	Stainless Steel wire Ø 8 mm Gallery height: 50 mm Gallery suitable for all shelves of a serving trolley	
0114503	For shelf-dimension: 800 x 500 mm Gallery 3-sided for SW/SSW-100 with shelf size 1000/600 Stainless Steel 50 mm total height, 40 mm above shelf	120,00 € / Pc
		
	Stainless Steel wire Ø 8 mm Gallery height: 50 mm Gallery suitable for all shelves of a serving trolley	
0114507	For shelf-dimension: 1000 x 600 mm Gallery 4-sided for SW/SSW-60 with shelf size 600/400 Stainless Steel 50 mm total height, 40 mm above shelf	129,00 € / Pc
		

Stainless Steel wire Ø 8 mm  
 Gallery height: 50 mm  
 Gallery suitable for all shelves of a serving trolley

- 0114504 For shelf-dimension: 600 x 400 mm  
 Gallery 4-sided for SW/SSW-80 with shelf size 800/500 Stainless Steel 50 mm total height, 40 mm above shelf 129,00 € / Pc



Stainless Steel wire Ø 8 mm  
 Gallery height: 50 mm  
 Gallery suitable for all shelves of a serving trolley

- 0114505 For shelf-dimension: 800 x 500 mm  
 Gallery 4-sided for SW/SSW-100 with shelf size 1000/600 Stainless Steel 50 mm total height, 40 mm above shelf 129,00 € / Pc



Stainless Steel wire Ø 8 mm  
 Gallery height: 50 mm  
 Gallery suitable for all shelves of a serving trolley

- 0116192 For shelf-dimension: 1000 x 600 mm  
 Waste Container Stainless Steel with hanging brackets 420 x 200 x 250 mm 155,00 € / Pc



- 0114802 Suitable for 8 x 5 and 10 x 6  
 Waste Bag Holding Frame - MHR 120 - For SW 8 x 5 and SW 10 x 6 176,00 € / Pc



To keep waste bags 120 litres.

- 0114584 External dimensions: 350 x 283 x 736 mm  
 Reinforced bar for SW/SSW-60 with shelf size 600/400 Stainless Steel 188,00 € / Pc



From high-quality Stainless Steel  
 Shelf load increase = 20 kg per shelf.

- 0114582 Reinforced bar for SW/SSW-80 with shelf size 800/500 Stainless Steel 188,00 € / Pc



From high-quality Stainless Steel  
Shelf load increase = 20 kg per shelf.

0114583 Reinforced bar for SW/SSW-100 with shelf size 1000/600 Stainless Steel 188,00 € / Pc



From high-quality Stainless Steel  
Shelf load increase = 20 kg per shelf.

### Transport Trolleys -> Basket Trolleys, Platform Trolleys, Transport Trolleys, Single Ended Trolleys -> Articles

Order No. article €

7500196 Basket trolley KR/ 60-60 770/675/855 Stainless Steel 460,00 € / Pc



Dimensions W x D x H mm:	770 x 675 x 855 mm
-----------------------------	--------------------

For Sterile goods baskets according to DIN  
Trolley manufactured from high-quality  
Stainless Steel

Stable, mobile round tube construction

Mobile on 4 swivel castors, 2 of them with brake,  
in galvanised-chromed housings, Ø 125 mm, with  
double ball-bearing, electro-conductive.  
2 bumpers protect trolley against damage.

7504872 Basket dispenser OKO/ 61-61 StU Stainless Steel to take 2 pcs. StU-baskets side by side 1500,00 € / Pc



Dimensions W x D x H mm:	715x830x900
-----------------------------	-------------

7504874 Basket Dispenser OKO / 61-41 ISO/SPRI to take ISO Norm baskets Stainless steel 1350,00 € / Pc



Dimensions W x D x H mm:	715x805x900
-----------------------------	-------------

0112602 Platform Trolley PW/8x5 Stainless Steel 498,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Deep-drawn metal platform with elevated profile edge all around
Platform-dimensions:	800 x 500 mm

Bumping-protection:	4 bumpers mounted above the castors made from elastic and shock-resistant plastic
Chassis:	4 swivel castors, 2 of them with total brake, Ø 125 mm
Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	120 kg
Net weight:	20 kg
External dimensions W x D x H in mm:	900 x 600 x 1020 mm

0112632 Platform Trolley PW/10x6 Stainless Steel

541,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Deep-drawn metal platform with elevated profile edge all around
Platform-dimensions:	1000 x 600 mm
Bumping-protection:	4 bumpers mounted above the castors made from elastic and shock-resistant plastic
Chassis:	4 swivel castors, 2 of them with total brake, Ø 125 mm
Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	120 kg
Net weight:	24 kg
External dimensions W x D x H in mm:	1100 x 700 x 1020 mm

0112662 Platform Trolley PW/5,4x5,4 Stainless Steel

526,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Flat metal platform, folded down by 80 mm all around
Platform-dimensions:	545 x 545 mm
Bumping-protection:	2 plastic bumpers at the front
Chassis:	4 swivel castors, 2 of them with total brake, Ø 125 mm
Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	120 kg
Net weight:	16 kg
External dimensions W x D x H in mm:	740 x 600 x 1020 mm

0112692 Platform Trolley PW/8,2x5,4 Stainless Steel

535,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Flat metal platform, folded down by 80 mm all around
Platform-dimensions:	820 x 545 mm
Bumping-protection:	2 plastic bumpers at the front
Chassis:	4 swivel castors, 2 of them with total brake, Ø 125 mm

Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	120 kg
Net weight:	20 kg
External dimensions W x D x H in mm:	1015 x 600 x 1020 mm

0162578 Heavy Duty Platform Trolley TWs/12x6 Stainless Steel

890,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Flat metal platform, folded down by 60 mm all around, additionally with reinforcing bar
Platform-dimensions:	1200 x 600 mm
Bumping-protection:	2 plastic bumpers at the front
Chassis:	2 swivel castors with total brake, 2 fixed castors, Ø 160 mm
Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	500 kg
Net weight:	28 kg
External dimensions W x D x H in mm:	1370 x 660 x 1050 mm

0162588 Front panelling for TWs/12 x 6 Stainless Steel

86,00 € / Pc



Front panelling manufactured from high-quality Stainless Steel, 1,0 mm thick, edged by 21 mm, rolled to the inside.

0114192 Dimensions: 552 x 526 mm Single Ended Trolley STWP/9,8x5 with welded cover plate

1006,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Platform manufactured from round tube, Ø 25 mm, with welded-on stainless steel plate
Platform-dimensions:	850 x 500 mm
Chassis:	2 swivel castors with total brake, 2 fixed castors, Ø 200 mm
Castor specification:	Corrosion-protected and maintenance-free

Load capacity:	500 kg
Net weight:	42 kg
External dimensions W x D x H in mm:	970 x 500 x 940 mm

0114212 Single Ended Trolley STWP/11,3x7 with welded cover plate

1159,00 € / Pc



Material:	High-quality Stainless Steel
Platform-type:	Platform manufactured from round tube, Ø 25 mm, with welded-on stainless steel plate
Platform-dimensions:	1000 x 700 mm
Chassis:	2 swivel castors with total brake, 2 fixed castors, Ø 200 mm
Castor specification:	Corrosion-protected and maintenance-free
Load capacity:	500 kg
Net weight:	52 kg
External dimensions W x D x H in mm:	1120 x 700 x 940 mm

General Conditions of Sale and Delivery / General Conditions