


IQAir CleanZone H13

The IQAir® CleanZone H13 is a table-top workstation



The IQAir® CleanZone H13

is a table-top workstation with double-function: clean air hood and particle extraction hood system for the professional high-efficiency filtration of bacteria, viruses, spores, dust, soot, smoke particles, allergens, aerosols and other particulate contaminants.

Applications

Cleanroom facilities, laboratories, medical & pharmaceutical environments, military installations, biohazard control units, government institutions, schools, universities, mail rooms.

All systems come complete with filters, owners manual, set of 4 casters, 10 foot power cord, and certificate of guaranteed performance.

Standard Special Features

Ultra-Quiet Design
Easy Filter Replacement
No-Draft Design
Long Warranty

Advanced Controls
Long Filter Life
Energy Efficient

Certified Performance
Triple Seal Design
High-Quality Fan
Tough Construction

Specifications

Performance Characteristics

Filtration-efficiency at 0.3 microns: greater than 99.97% (certified)
in m3/h (Speed 1-5): 60, 110, 160, 240, 380 (certified)

Material/Colour:

Hood: acrylic (transparent)
Housing: steel (light blue, powder-coated) & ABS (light grey)

System includes filters, hood, power cable, instruction manual and certificate of performance.

Dimensions/Weight:

Total (DxWxH): 56 x 72 x 90 cm / 18,5 kg (able to be taken apart)
Housing: 42 x 37 x 45 cm / 12,8 kg
Hood: 56 x 72 x 47 cm / 5,7 kg

CR H13 HEPA-Filter (individually tested & certified)

Type: cleanroom grade HEPA filter;
Class: H13 (L) Efficiency: greater than 99.97% at 0.3 microns

Surface area: ca. 4.1 m²

PreMax™ Filter

Type: fine-dust pre-filter;
Class: F8 (S) Efficiency: greater than 55% at 0.3 µm
Surface area: ca. 2.3 m²


Other Information

For more filtration of microbiological contaminants, bacteria, and viruses see our [Cleanroom Series](#) or [Dental Series](#) models.

For more critical control of gaseous compounds see our [GC Series](#) models.



IQAir CleanZone H13
shown in a mail room environment



Particle Extraction Hood



Clean Air Hood



HOSPITAL ENGINEERING CZ

Člen Asociace výrobců a dodavatelů zdravotnických prostředků

HOSPITAL ENGINEERING CZ, s.r.o.
Nádražní 7, 602 00 Brno
tel./fax.: 00 420 543 331 101-2
www.hospitalengineering.cz