



Colon Hydromat comfort: for hydro-colon therapy

Colon Hydromat comfort has been equipped with an integrated colour control. Neatly arranged control elements and a clear design guarantee for simple handling.



Highlights of this system:

- More ease of use and more security thanks to automatic ventilation
- An integrated water flow counter allows for a better overview of the process
- Multiple security mechanisms guarantee for highest safety
- Precise notification about treatment status on a display
- Integrated colour control

Technical data

Name:	Colon Hydromat comfort
Power supply:	100-240 V~50-60 Hz
Water pressure (cold and hot):	2-4 bar
Size (WxHxD):	510 x 545 x 250 mm
Weight:	27.5 kg
Classification according to Directive 93/42 EEC	II A

- Ord. no. 15701 left
- Ord. no. 15702 right
- Ord. no. 15703 left
- Ord. no. 15704 right