


**Treatment options with oxygen/ozone mixture in their latest form.**  
**Every user will be delighted at the new design as well as the user's guidance via touchscreen.**

Suitable for:  
the hyperbaric and the normobaric therapy  
ozone taking by syringe permanent gassing  
intestinal insufflation



Ozone trolley FUEGO optional available

## Highlights of this system:

- easy to use device with touchscreen and integrated user's guidance
- Input of concentration with number keys infinitely adjustable
- The entered concentration is automatically measured and controlled by an analyser
- All connections are Luer-Lock

## Technical data

Name of appliance:	Hyper Medozon comfort
Power supply:	230V~/50 Hz oder 115V~/50 Hz
Power input:	310 VA
Ozone concentration:	0 und 5-80 µg/ml (control via analyser)
Flow rate:	1,0 l/min
Dimensions (BxHxT):	535 x 180 x 460 mm
Weight:	18,7 kg ( with airtrap )
Product category acc. to MP-directives 93/42 EWG:	II b
	■ 230V Ord. no. 16640
	■ 115V Ord. no. 16641
	■ 230V Ord. no. 16643
	■ 115V Ord. no. 16644