



Hyper Vetozone comfort: normobaric and hyperbaric ozone therapy for pets

Hyper Vetozone comfort enthralls with its intuitive handling. The system can be used for normobaric as well as for hyperbaric ozone therapies. The touch screen control allows adjusting the concentration of the ozone gas quickly and comfortably.



Highlights of this system:

- Can be used in normobaric and hyperbaric treatments
- High flexibility thanks to the possibility to draw up ozone into a syringe
- Implementation of permanent external ozone treatment
- Multilingual operational guidance via touch screen in colour
- Reinfusion pressure of up to 1.2 bar
- High compatibility thanks to Luer Lock connections

Technical data

Name:	Hyper Vetozone comfort
Power supply:	230V~/50 Hz oder 115V~/50 Hz
Input:	210 VA
Concentration ozone:	5-80 µg/ml
Flow through:	1 l/min
Size (BxHxT):	535 x 180 x 460 mm
Weight:	18.7 kg

■ Ord. no. 25000

Hyper Vetozone comfort