

Magnetic field therapy



Trimagx: system for magnetic field therapy

- Pulsed magnetic waves
- High flexibility thanks to annular coils in three different sizes
- Touch screen control and a lucid design increase clarity
- Fully automatic treatment facilitates operational guidance
- Its compact size qualifies the systems for use in every doctor's surgery
- Multilingual operational guidance



Technical data

Name:	Trimagx
Power supply:	230V~/50 Hz
Input:	470 mA
Size (WxHxD):	535 x 1130 x 770 mm
Size of coils:	550 mm, 350 mm, 200 mm
Weight:	78 kg
Classification according to Directive 93/42 EEC	II a

■ Ord. no. 14002

Trimagx