



Vetozone compact: for the ozone therapy in the veterinary sector

Excellent operational comfort and high safety – these are the characteristics of systems by HERRMANN Apparatebau. Vetozone compact has been developed according to our strict quality guidelines and can be used for almost the entire range of modern ozone treatment. The integrated colour touch screen allows for an uncomplicated handling.



Highlights of this system:

- Can be used for all types of normobaric ozone therapy
- High flexibility thanks to the possibility to draw up ozone into a syringe
- High compatibility thanks to Luer Lock connections
- Multilingual operational guidance (German, English, Spanish, Italian, Turkish)
- Colour touch screen for a more comfortable handling

Technical data

| | | | |
|----------------------|----------------------------|--------|-----------------|
| Devicename: | Vetozone compact | | |
| Power supply: | 230V~/50 Hz or 115V~/50 Hz | | |
| Input: | 240 VA | | |
| Concentration ozone: | 0;5-80 µg/ml | | |
| Flow through: | 1,0 l/min | | |
| Dimensions (WxHxD): | 424 x 185 x 388 mm | ■ 230V | Order no. 25002 |
| Weight: | 13,1 kg | ■ 115V | Order no. 25003 |